Neu ab:
New since:
Nouveau depuis:
Nuevo desde:


01.2004

sirona
The Dental Company

ORTHOPHOS Plus DS/Plus DS Ceph

LAN Kit A: 59 00 357
LAN Kit B: 59 17 401
LAN Kit C: 58 98 288

Einbauanleitung für Nachrüstsatz LAN-Anbindung
Installation Instructions for LAN interface upgrade kit
Instr. de montage pour kit d'extension "connexion LAN"
Instr. de montaje para el kit de ampliación de conexión LAN

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
sirona
The Dental Company

ORTHOPHOS Plus DS/Plus DS Ceph

Installation instructions
for LAN interface upgrade kit

LAN Kit A: 59 00 357
LAN Kit B: 59 17 401
LAN Kit C: 58 98 288

Valid for units
from Serial No.: 1,000/50,000

Contents

1 Safety....5

2 General notes....6

Part A

3 LAN Kit A, ORTHOPHOS Plus DS 96/98 / 59 00 357....7

4 (A) Parts supplied....8

5 (A) Preparatory work....9

6 (A) Dismantling the old electronics carrier 10

7 (A) Installing the new electronics carrier...14

8 (A) Connecting additional protective ground wires....17

9 (A) Installing and configuring Orthophos Plus DS....18

10 (A) Adjusting board DX1....19

11 (A) Generating a test image (Panoramic unit)....20

12 (A) Generating a test image (Ceph unit)...22

13 (A) Phantom exposure....23

Part B

14 (A) Final Work....24

15 LAN Kit B, ORTHOPHOS Plus DS (+DX1) / 59 17 401....25

16 (B) Parts supplied....26

17 (B) Preparatory work....27

18 (B) Replacing board DX1....28

19 (B) Installing the LAN slide-in module....29

20 (B) Connecting additional protective ground wires....30

english

21 (B) Installing and configuring Orthophos Plus DS 31

22 (B) Adjusting board DX1 32

23 (B) Generating a test image (Panoramic unit) 33

Part C

24 (B) Generating a test image (Ceph unit) ... 35

25 (B) Phantom exposure 36

26 (B) Final work 37

27 LAN Kit ORTHOPHOS Plus DS / 58 98 288 38

28 (C) Parts supplied 39

29 (C) Preparatory work 40

30 (C) Installing the LAN slide-in module 41

31 (C) Connecting additional protective ground wires 42

32 (C) Installing and configuring Orthophos Plus DS 43

33 (C) Generating a test image (Panoramic unit) 44

34 (C) Generating a test image (Ceph unit) ... 46

35 (C) Phantom exposure 47

36 (C) Final Work 48

1 Safety

ATTENTION

It is essential that you comply with the warning and safety notes contained in the Installation Instructions.

They are specifically identified by NOTE, ATTENTION or WARNING.

ATTENTION

The installation must be carried out in accordance with the requirements stated in our Installation Instructions.

ATTENTION

The installation may be carried out only by personnel specifically authorized by Sirona.

ATTENTION

It is essential that you switch off the device PRIOR to beginning the installation.

ATTENTION

Any person who assembles or modifies a medical electrical system complying with the standard EN 60 601-1-1 (Safety requirements for electromedical systems, Section 2.203) by combination with other equipment is responsible for ensuring that the requirements of this regulation are fulfilled to their full extent for the safety of the patients, the operators and the environment.

ATTENTION


To guarantee the operational safety of electromedical devices, the use of mobile telephones is prohibited in the practice or hospital area.

ATTENTION

When opening the unit:
Please observe precautionary measures when handling PC boards (ESD).
Touch a ground point to discharge static electricity before handling any components.

NOTE

Cable ties
When cutting through cable ties with side cutters make sure that no cables are cut through!

english

2 General notes

These installation instructions describe the conversion of an ORTHOPHOS Plus DS / Plus DS Ceph to a network-capable unit (LAN).

Depending on the serial number of the unit to be converted, the work which needs to be carried out is different.

The installation instructions are therefore divided into three parts: **A**, **B** and **C**.

For ease of use, the relevant letter is placed in front of the chapter headings.

Part	Serial No. Panoramic unit	Serial No. Ceph unit	Conversion action	Part No.	Page
A	1,000 – 1,399	50,000 – 50,399	Replacing the electronics carrier	59 00 357	7 - 24
B	1,400 – 1,445	50,400 – 50,459	Replacing the board DX1 and replacing the board DEB by the LAN slide-in module	59 17 401	25 - 37
C	1,446 – 1,999	50,460 – 51,999	Replacing the board DEB by the LAN slide-in module	58 98 288	38 - 48

NOTE

Sidexis Version 5.2 or higher as well as the ORTHOPHOS Plus DS Service Manual is required for the installation.

NOTE

For the LAN Kit A (replacing the electronics carrier) you require apart from the customary tools a ball socket wrench (SW 5), in order to loosen or tighten the mounting screws of the electronics carrier which are difficult to access.

NOTE

After conversion, the unit is operated through a network using a LAN cable. This requires appropriate planning work in advance in accordance with
"SIDEXIS, Installation Planning" ,
Order Number 58 22 528,
and preliminary work must be carried out in accordance with
"SIDEXIS, Installation Stand-alone Systems",
Order Number 58 72 952, or
"SIDEXIS, Software Service Manual, Installation of Multistation Systems“ *,*
Order Number 58 73 208
.

NOTE

*The **XOP** plug-in board, which is not required in the PC after the conversion, can be removed. It can be returned to the factory for credit.*

ATTENTION


*If you switch the unit off after installing the LAN slide-in module, you must comply with a **minimum waiting time of 20 sec.** before you switch back on.*
Otherwise there is a risk that the XAB OP board contained in the LAN slide-in module is not initialized.

3 LAN Kit A, ORTHOPHOS Plus DS 96/98 / 59 00 357


This part of the instructions applies for the following unit versions:

- Panoramic units with serial numbers from **1,000 to 1,399**
- Ceph units with serial numbers from **50,000 to 50,399**

Part A

4 (A) Parts supplied


1. Technical documents
2. Nameplate
3. Membrane keyboard with mounting plate and mounting hardware (required only for older units)
4. Electronics carrier, pre-assembled
5. Grounding straps with mounting hardware (required only for older units; long: for ring carrier; short: for motor carrier)
6. Cable ties
7. Program memory card
8. Cap clip, screws (LAN cable strain relief)
9. LAN cable
10. Covers
11. Insulating foil for display board DX7
12. Cable protection 2x

Ceph unit

13. Take over the cables:
 L24 from DX1 X15 to DX91 X6 and
 L25 from DX1 X30 to DX 91 X4

5 (A) Preparatory work

NOTE

Note the user-specific settings of the unit (refer to Operating Instructions).

1. Switch the unit off and disconnect it from the power supply.


ATTENTION


Danger of electric shock. If the unit is opened too soon: The capacitors in the unit need at least 10 seconds to discharge.

2. Unscrew the multitimer or the remote control.
3. Remove the cable connection to the PC.
4. Loosen the four top screws on the front cover and release it.
5. Remove the lead cover (5 screws).
6. Remove the side cover (2 screws).


6 (A) Dismantling the old electronics carrier


Required only for older units with membrane keyboard without grounding strap!

Dismantling the door

1. Unscrew the cover plate of the power supply unit (2 screws).
2. Loosen the two screws of the plate of the ground connection (do not unscrew) and withdraw the grounding strap.
3. Remove the 6 screws **A** of the inner cover and unscrew the PC board located below.
4. Remove the upper 2 Allen screws of the door cover and take it off.
5. In older units with **membrane keyboard without grounding strap** this must be replaced:
 Loosen the old membrane keyboard carefully with a knife at the lower right corner and pull it off from the door. Affix the new membrane keyboard.

(A)

6. Pull the connector **X16** off from the board **DX1** and draw the cable back through the slot in the electronics plate. Unscrew the correction switch.
Re-installation!
The cable and correction switch must be installed again after the conversion.

7. Loosen the grounding strap for the motor carrier on the right side.


8. Now loosen all cables connected with the electronics carrier (cut the cable ties through if necessary):
DX1: X4
DX1: X5
DX1: X7
DX1: X9
DX1: X11
DEB: X6
DEB: X7

only for Ceph:
DEB: X3 / X4 (loosen lateral clamp)
DX91: X1 / X2 / X3 / X7 (loosen lateral clamp)


9. Loosen the plug connection DX31: X2 and draw the cable back (the cable is run through the electronics carrier).

10. Remove the board DX91 (only for Ceph).


(A)

Left


Right (DEB)

Front (DX1)

11. Remove the mounting screws of the column cover (apron)
- on the left side under the fuses,
- on the right side through the board DEB
- and 2 screws at the front through the board DX1.

(A)

Unscrewing the electronics carrier

The front mounting screws are difficult to access. You need a ball socket wrench (SW 5) to loosen the screws.

ATTENTION

Tape over the painted panel of the bite block with adhesive tape. There is a risk of scratching this when removing the electronics carrier!

12. Remove firstly the two lateral mounting screws of the electronics carrier (left-side transformer mounting).

13. Then remove the front mounting screws in the order shown
a, **b** and **c**


ATTENTION

Hold the electronics carrier firmly (it has sharp edges!) when you remove the last screw. Otherwise there is a risk that the carrier drops down.

14. Let the electronics carrier down carefully and pull it slightly obliquely inclined out to the front.

15. Take over the mounting bolts (2x left, 1x right), the rubber spacer from the old to the new electronics carrier.

7 (A) Installing the new electronics carrier


1. Fasten the new electronics carrier firstly with the front screws. Then screw in the side screws on the transformer (see page 13).
2. Align the apron and fasten on the left, right and front (see page 12).
 Take care that it does not rub against the column.
3. Reattach the door and screw it on tightly using the 2 screws.
4. Replace the display board DX7 in the door (already cabled with the electronics carrier). For this purpose fix the 4 red spacer O-rings of the old DX7 with a little grease to the inside of the door and fasten the insulating foil on the conductor side of the board.
5. Connect the grounding strap of the membrane keyboard.
 If the membrane keyboard was replaced previously, corresponding mounting holes must be drilled. Use the installation plate as drilling template.
6. Screw the door cover back on again.
 (6 screws).
7. Restore the ground connection to the motor carrier again (see page 11).


In older units the following grounding straps must be installed:

8. Grounding strap for the ring carrier (long strap):
 a) Remove the mirror (it is plugged on the hinge pin; grasp with both hands!)
 b) Unscrew the cover (2 screws)
 c) Fasten the grounding strap
 d) Fasten the cover again
 e) Plug the mirror back on
9. Grounding strap for the motor carrier (short strap):
 a) Unscrew the capacitor and screw it back on together with the grounding strap.
 b) Fasten the other end of the grounding strap to the inside of the carrier.

58 97 801 D 3297
D 3297.092.02.04.02

(A)

10.

11.


10. Insert the cable protection and restore all cable connections to the electronics carrier:
DX1: X4
DX1: X5
DX1: X7
DX1: X8
DX1: X9
DX1: X11
XAB OP: X6
XAB OP: X7

11. Only for Ceph:
Refit the board DX91 (see page 11)
and make the cable connections:
XAB OP: X3 / X4
(Fasten cable to upper clamp **A**)
DX91: X1 / X2 / X3 / X7
(Fasten cable to lower clamp **B**)
Take over the cables (see note on page 8):
L24 from DX1 X15 to DX91 X6 and
L25 from DX1 X30 to DX 91 X4.

NOTE
Fasten all required cable ties only after a successful test.

english

(A)

12. Install the cable to DX31: X2 (L10) again (through the electronics carrier).


13. Refit the correction switch (installation and adjustment see Service Manual) and run the cable through the slot in the electronics plate.
Plug **X16** back on board **DX1** (see page 11).

14. Check or set the line voltage by changing over the plug on the power supply unit board.

15. Insert the program memory card and secure it with the locking hook **R**.

16. Paste over the label on the left side of the column with the new software version.

8 (A) Connecting additional protective ground wires

NOTE


An additional protective ground wire (4mm² with cable lug 4-6 DIN 46234 CU) is required for the X-ray unit and the connected PCs (provided they stand in the same room). It must comply with the requirements for equipotential bonding.

ATTENTION

The electrical installation must be made according to VDE 0107 (VDE 0170, 10.94, Equipotential bonding: Section 4.4.1 Application group 1).

1. Connect the additional protective ground wire for the ORTHOPHOS Plus in the supply section on the floor plate in front of the supply terminals.
 Run the cable in the rear section through the sleeve **T** to the outside.
2. Connect the protective ground wire to the building ground. It is pre-assembled with a cable lug 5-2.5 DIN 46234. The cable lug can be removed if the cable is connected to a terminal.

9 (A) Installing and configuring Orthophos Plus DS

NOTE

SIDEXIS 5.3 or higher and the SIXABCon configuration program must be installed for the following work.


1. Run the SIXABCon configuration program, see SIDEXIS installation of single/multistation systems. (Otherwise the service routines S.03, S.04 and S.05 cannot be performed.)
2. Install the Pan/Ceph sensor floppy using the SIXABCon configuration program.
3. Set SIDEXIS to factory service.

ATTENTION

IMPORTANT NOTES

- *For the following tests, it is important to ensure compliance with the **radiation protection guidelines** (see Operating Instructions).*
- *The **radiation** remains on as long as the release button is pressed.*
- *As long as radiation is active the **X-RAY** radiation on/off indicator must be lit up.*
- *The radiation on-time is also indicated **by an acoustic signal** emitted by the unit.*
- *The prescribed **cooling periods** are monitored by the **automatic exposure blocking** (see Operating Instructions).*

10 (A) Adjusting board DX1

- Reconnect the unit to the power supply and switch it on.

After switching ON the self-adjustment runs in the unit. Wait until the customary displays for program/exposure time, kV and mA light up on the multitimer.

Make the following checks and settings according to the Service Manual.

1. Panoramic unit.

- Service routine **S.03**:
 Set **nominal values** for kV, mA and preheating.
 Perform test step **04** completely.
- Service routine **S.05**:
 Select test step **02** (automatic).
 Perform **heating adaptation**.
- Select service routine **S.04** test step **03**.
 Check **actual values** for kV, mA and preheating.
- Select service routine **S.21**.
 Program **diaphragms**.
- Select service routine **S.20**.
 Set **forehead support width**.
- Select service routine **S.18**.
 Check **height adjustment** for function (PWM).
- Select service routine **S.25**.
 Set kV/mA level series according to national regulations and possibly customer's wishes.
- To leave the service routines switch the unit OFF and back ON.


2. Ceph unit.

- Select service routine **S.34** test step **01** and **03**.
 Check the Ceph arm.
- Select service routine **S.18** test step **07**.
 Adapt the PWM, speed adaptation for Ceph movement.

The board DX1 is now set correctly for the X-ray tube assembly.

english

11 (A) Generating a test image (Panoramic unit)

NOTE

Do not press the test button on the multitimer before generating the test image.

1. Start SIDEXIS.
2. Activate the **Service** menu.
3. Activate the **Constancy test** menu item.
4. In the **X-ray components** menu, activate the **New** menu window.
5. Enter the service password.
 (see SIDEXIS Service Manual).

The **New X-ray component** window will appear.

6. Fill in all the input fields.
7. Press the Return key to continue. If you do not know all the information, enter **Unknown** or **?** .
8. Press the **OK** button.
9. Activate the **XC XP exposure performance** icon**.**

This displays the **Select test type** selection dialog:


- Service exposure
- Acceptance test
- Constancy test
- Cancel

10. If another dialog is displayed:
 Cancel
 and activate **XC XP**.
11. From the **Select test type** dialog, select
 Service exposure.

The selection dialog **Select service exposure** appears.

12. From this dialog, select **Digital test pattern (1)**.

(A)

13. Press the release button.

This triggers a cycle **with no** release of radiation.

A dialog box then appears.

If **OK** appears on the display, the data path is in order.

14. Press the **OK** button.

The test image generated is displayed on the screen.


If **Error** appears on the display, please refer to the Service Manual.

15. Close the test image ([Ctrl]+[F4]).

16. Press the Return key **R**.

The X-ray tube unit moves to its start position.

12 (A) Generating a test image (Ceph unit)

Checking the position of the patient diaphragm:

1. Set the primary diaphragm **30** on the diaphragm wheel (locking button **D**).
2. Press the Return key **R**.
 The X-ray tube assembly moves into the position for cephalography.
 The patient diaphragm moves into the position for asymmetric exposures.
3. Set the primary diaphragm **40** on the diaphragm wheel (locking button **D**).
4. Press the Return key **R**.
 The patient diaphragm moves into the position for symmetric exposures.

- If necessary adjust the diaphragms according to the installation instructions.

Generate test image:

- From the
 Select service exposure dialog activate the option **Digital test pattern (1)**.

5. Press the release button.
 This triggers a scan with no release of radiation.

A dialog box then appears.
If

OK

appears on the display, the data path is OK.

Press the **OK** button.
The test image generated is displayed on the screen.


If

Error

appears on the display, see the information in the Service Manual.

- Close the test image.

13 (A) Phantom exposure

1. Insert the needle phantom as far as it will go.

2. Set the **primary diaphragm 10** on the diaphragm wheel.

3. Press the Return key **R**.

The X-ray tube unit moves to its start position.

4. From the selection dialog box on the screen, select **Service exposure** and activate **Quality test exposure (4)**.


5. Release an exposure for a **complete cycle** of the rotating part.

6. Evaluate the image on the screen.
Measure the line spacings using **Measure length** in the **Analyze** menu (see SIDEXIS Operator's Manual, Functional Description, Image processing, Measure length).

- If the tolerance given in the diagram is exceeded, you must adjust the M2/M3 actuators (see Service Manual).
- Exit exposure ([Ctrl]+[F4]).

7. Take the needle phantom back out again.

14 (A) Final Work

1. Perform the national acceptance test.
2. Reattach all separated cable ties. Secure the LAN cable with the strain relief.
3. Affix the nameplate to the right side of the apron.
4. Screw the lead cover back on again.

ATTENTION

*For reasons of electromagnetic compatibility, **all** screws must be turned down and tightened.*

5. Attach the new side cover using the two screws.
6. Attach the new front cover and screw it on tightly using the four screws.
7. Screw the cover of the power supply unit back on.


- Restore the previously noted user-specific settings of the unit (refer to Operating Instructions).

15 LAN Kit B, ORTHOPHOS Plus DS (+DX1) / 59 17 401


This part of the instructions applies for the following unit versions:

- Panoramic units with serial numbers from **1,400 to 1,445**
- Ceph units with serial numbers from **50,400 to 50,459**

Part B


16 (B) Parts supplied

1. Technical documents
2. LAN slide-in module
3. Program memory card, V035 or higher
4. Cable ties
5. Cables
 a) LAN
 b) L32
 c) L21
 d) Protective ground wire
6. Board DX1
7. Board DX32

17 (B) Preparatory work

NOTE

Note the user-specific settings of the unit (refer to Operating Instructions).


1. Switch the unit off and disconnect it from the power supply.

ATTENTION

Danger of electric shock. If the unit is opened too soon: The capacitors in the unit need at least 10 seconds to discharge.


2. Disconnect the multitimer.
3. Remove the cable connection to the PC.
4. Loosen the four top screws on the front cover and release it.
5. Remove the lead cover (5 screws).
6. Remove the side cover (2 screws).
7. Remove the connector of cable L30, the cover plate containing the recess for the PC connection on the DEB board (no longer required) on the back of the unit. You will need to keep the screws to attach the XAB slide-in module.

18 (B) Replacing board DX1

1. Now loosen all cables connected with the board **DX1**
 (cut the cable ties through if necessary):
 DX1: X1 (at DX32:X3; L4)
 DX1: X2 (at DX32:X4; L5)
 DX1: X4 (at K3 / K9; L10)
 DX1: X5 (at K3; L10)
 DX1: X7 (at DX5: X1; L13)
 DX1: X8 (at DX7: X2; L6)
 DX1: X9 (at K7 / K8; L7)
 DX1: X10 (at K2; L3)
 DX1: X11 (at DX5: X2; L15)
 DX1: X14 (at DEB: X14; L17)
 DX1: X15 (at DX91: X6; L24; only Ceph)
 DX1: X16 (at correction switch S_{HV})
 DX1: X30 (at DX91: X4; L25; only Ceph)
2. Unscrew the board (5 plastic screws) and replace it by the new board.

19 (B) Installing the LAN slide-in module

ATTENTION

While installing the PC boards use the grounding arm band present in the unit.


1. Loosen all plug connections on the **DEB** board: X5, X6, X7, X14 (Ribbon cable), for ORTHOPHOS Plus DS Ceph in addition X3 and X4.
2. Unscrew the **DEB** board (4 plastic nuts with spacer sleeves and plastic plate - no longer required). Remove backing nut from DEB grounding stud.
3. Fasten the LAN slide-in module using the screws from the cover plate (see page 40 section 7.).
4. Replace the cable L19 DS with the cable L32. Connect cable L32 to board **XAB OP X5** and to the power supply board **DX32 X7**; if no X7 is present on the installed DX32, install the supplied new DX32.
5. Connect power supply board **DX32: X6** with the board **XAB OP: X10**. Use the cable L21 / No. 33, DS for this purpose.
6. Plug the cables X6, X7, X14 on the board **XAB OP** (for ORTHOPHOS Plus DS Ceph in addition X3 and X4).
7. Restore all plug connections on the board **DX1** (see page 28).
8. Plug the LAN cable on the board **XAB-OP** and make the cable connection to the network.

NOTE

Fasten all required cable ties only after a successful test.

9. Insert the program memory card and secure it with the locking hook **R**.
10. Paste over the label on the left side of the column with the new software version.

20 (B) Connecting additional protective ground wires

NOTE


An additional protective ground wire (4 mm² with cable lug 4-6 DIN 46234 CU) is required for the X-ray unit and the connected PCs (provided they stand in the same room). It must comply with the requirements for equipotential bonding.

ATTENTION

The electrical installation must be made according to VDE 0107 (VDE 0170, 10.94, Equipotential bonding: Section 4.4.1 Application group 1).

1. Connect the additional protective ground wire for the ORTHOPHOS Plus in the supply section on the floor plate in front of the supply terminals.
 Run the cable in the rear section through the sleeve **T** to the outside.
2. Connect the protective ground wire to the building ground. It is pre-assembled with a cable lug 5-2.5 DIN 46234. The cable lug can be removed if the cable is connected to a terminal.
3. Secure the protective ground wire with a cable tie as strain relief to the supply cable.

21 (B) Installing and configuring Orthophos Plus DS

NOTE

SIDEXIS 5.3 or higher and the SIXABCon configuration program must be installed for the following work.


1. Run the SIXABCon configuration program, see SIDEXIS installation of single/multistation systems. (Otherwise the service routines S.03, S.04 and S.05 cannot be performed.)
2. Install the Pan/Ceph sensor floppy using the SIXABCon configuration program.
3. Set SIDEXIS to factory service.

ATTENTION

IMPORTANT NOTES

- *For the following tests, it is important to ensure compliance with the **radiation protection guidelines** (see Operating Instructions).*
- *The **radiation** remains on as long as the release button is pressed.*
- *As long as radiation is active the **X-RAY** radiation on/off indicator must be lit up.*
- *The radiation on-time is also indicated **by an acoustic signal** emitted by the unit.*
- *The prescribed **cooling periods** are monitored by the **automatic exposure blocking** (see Operating Instructions).*

22 (B) Adjusting board DX1

- Reconnect the unit to the power supply and switch it on.

After switching ON the self-adjustment runs in the unit. Wait until the customary displays for program/exposure time, kV and mA light up on the multitimer.

Make the following checks and settings according to the Service Manual.

1. Panoramic unit.


- Service routine **S.03**:
 Set **nominal values** for kV, mA and preheating.
 Perform test step **04** completely.
- Service routine **S.05**:
 Select test step **02** (automatic).
 Perform **heating adaptation**.
- Select service routine **S.04** test step **03**.
 Check **actual values** for kV, mA and preheating.
- Select service routine **S.21**.
 Program **diaphragms**.
- Select service routine **S.20**.
 Set **forehead support width**.
- Select service routine **S.18**.
 Check **height adjustment** for function (PWM).
- Select service routine **S.25**.
 Set kV/mA level series according to national regulations and possibly customer's wishes.
- To leave the service routines switch the unit OFF and back ON.

2. Ceph unit.

- Select service routine **S.34** test step **01** and **03**.
 Check the Ceph arm.
- Select service routine **S.18** test step **07**.
 Adapt the PWM, speed adaptation for Ceph movement.

The board DX1 is now set correctly for the X-ray tube assembly.

23 (B) Generating a test image (Panoramic unit)

NOTE

Do not press the test button on the multitimer before generating the test image.

1. Start SIDEXIS.
2. Activate the **Service** menu.
3. Activate the **Constancy test** menu item.
4. In the **X-ray components** menu, activate the **New** menu window.
5. Enter the service password.
 (see SIDEXIS Service Manual).

The **New X-ray component** window will appear.

6. Fill in all the input fields.
7. Press the Return key to continue. If you do not know all the information, enter **Unknown** or **?** .
8. Press the **OK** button.
9. Activate the **XC XP exposure performance** icon.

This displays the **Select test type** selection dialog:

- Service exposure
- Acceptance test
- Constancy test
- Cancel


10. If another dialog is displayed:
 Cancel
 and activate **XC XP**.
11. From the **Select test type** dialog, select
 Service exposure.

The selection dialog **Select service exposure** appears.

12. From this dialog, select **Digital test pattern (1)**.

(B)

13. Press the release button.

This triggers a cycle **with no** release of radiation.

A dialog box then appears.

If **OK** appears on the display, the data path is in order.

14. Press the **OK** button.

The test image generated is displayed on the screen.


If **Error** appears on the display, please refer to the Service Manual.

15. Close the test image ([Ctrl]+[F4]).

16. Press the Return key **R**.

The X-ray tube unit moves to its start position.

24 (B) Generating a test image (Ceph unit)

Checking the position of the patient diaphragm:

1. Set the primary diaphragm **30** on the diaphragm wheel (locking button **D**).
2. Press the Return key **R**.
 The X-ray tube assembly moves into the position for cephalography.
 The patient diaphragm moves into the position for asymmetric exposures.
3. Set the primary diaphragm **40** on the diaphragm wheel (locking button **D**).
4. Press the Return key **R**.
 The patient diaphragm moves into the position for symmetric exposures.

- If necessary adjust the diaphragms according to the installation instructions.

Generating a test image:

- From the **Select service exposure** dialog activate the option **Digital test pattern (1)**.

5. Press the release button.
 This triggers a scan with no release of radiation.

A dialog box then appears.
If

OK

appears on the display, the data path is OK.

Press the **OK** button.
The test image generated is displayed on the screen.


If

Error

appears on the display, see the information in the Service Manual.

- Close the test image.

25 (B) Phantom exposure

1. Insert the needle phantom as far as it will go.

2. Set the **primary diaphragm 10** on the diaphragm wheel.

3. Press the Return key **R**.

The X-ray tube unit moves to its start position.

4. From the selection dialog box on the screen, select **Service exposure** and activate **Quality test exposure (4)**.


5. Release an exposure for a **complete cycle** of the rotating part.

6. Evaluate the image on the screen.
Measure the line spacings using **Measure length** in the **Analyze** menu (see SIDEXIS Operator's Manual, Functional Description, Image processing, Measure length).

- If the tolerance given in the diagram is exceeded, you must adjust the M2/M3 actuators (see Service Manual).
- Exit exposure ([Ctrl]+[F4]).

7. Take the needle phantom back out again.

26 (B) Final work

1. Perform the national acceptance test.
2. Reattach all separated cable ties. Secure the LAN cable with the strain relief.
3. Affix the nameplate with the new software version to the right side of the apron.
4. Screw the lead cover back on again.

ATTENTION

*For reasons of electromagnetic compatibility, **all** screws must be turned down and tightened.*

NOTE

Check that the side cover (59 68 250 D3200) has a recess for the LAN connection. If this is not the case, you must saw out/file the recess according to the above sketch.

5. Attach the side cover again using the two screws.
6. Attach the new front cover and screw it on tightly using the four screws.

- Restore the previously noted user-specific settings of the unit (refer to Operating Instructions).

27 LAN Kit ORTHOPHOS Plus DS / 58 98 288


This part of the instructions applies for the following unit versions:

- Panoramic units with serial numbers from **1,446 to 1,999**
- Ceph units with serial numbers from **50,460 to 51,999**

Part C


28 (C) Parts supplied

1. Technical documents
2. LAN slide-in module
3. Program memory card, V035 or higher
4. Cable ties
5. Cables
 a) LAN
 b) L32
 c) L21
 d) Protective ground wire

29 (C) Preparatory work

1. Switch the unit off and disconnect it from the power supply.

ATTENTION

Danger of electric shock. If the unit is opened too soon: The capacitors in the unit need at least 10 seconds to discharge.

2. Disconnect the multitimer.

3. Remove the cable connection to the PC.


4. Loosen the four top screws on the front cover and release it.

5. Remove the lead cover (5 screws).

6. Remove the side cover (2 screws).

7. Remove the connector of cable L30, the cover plate containing the recess for the PC connection on the DEB board (no longer required) on the back of the unit. You will need to keep the screws to attach the LAN slide-in module.

30 (C) Installing the LAN slide-in module


ATTENTION

While installing the PC boards use the grounding arm band present in the unit.

1. Loosen all plug connections on the board **DEB**: X5, X6, X7, X14 (Ribbon cable), for ORTHOPHOS Plus DS Ceph in addition X3 and X4.
2. Unscrew the **DEB** board (4 plastic nuts with spacer sleeves and plastic plate - no longer required). Remove backing nut from DEB grounding stud.
3. Fasten the LAN slide-in module using the screws from the cover plate (see page 40 section 7.).
4. Replace the cable 19DS by the cable 32, DS (power supply unit board **DX32: X7**) and plug it on to X5 on board **XAB OP**.
5. Connect the power supply unit board **DX32: X6** with the board **XAB OP: X10**. Use the cable L21 / No. 33, DS for this purpose.
6. Plug the cables X6, X7, X14 on the board **XAB OP** (for ORTHOPHOS Plus DS Ceph in addition X3 and X4).
7. Plug the LAN cable on to the board **XAB-OP** and make the cable connection to the network.
 Secure the LAN cable with the strain relief.
8. Replace the program memory card and secure it with the locking hook **R**.
 Observe the notes in the instructions for **replacing the memory card** (54 52 854).
9. Paste over the label on the left side of the column with the new software version.

31 (C) Connecting additional protective ground wires

NOTE


An additional protective ground wire (4 mm² with cable lug 4-6 DIN 46234 CU) is required for the X-ray unit and the connected PCs (provided they stand in the same room). It must comply with the requirements for equipotential bonding.

ATTENTION

The electrical installation must be made according to VDE 0107 (VDE 0170, 10.94, Equipotential bonding: Section 4.4.1 Application group 1).

1. Connect the additional protective ground wire for the ORTHOPHOS Plus in the supply section on the floor plate in front of the supply terminals.
 Run the cable in the rear section through the sleeve **T** to the outside.
2. Connect the protective ground wire to the building ground. It is pre-assembled with a cable lug 5-2.5 DIN 46234. The cable lug can be removed if the cable is connected to a terminal.
3. Secure the protective ground wire with a cable tie as strain relief to the supply cable.

32 (C) Installing and configuring Orthophos Plus DS

NOTE

SIDEXIS 5.3 or higher and the SIXABCon configuration program must be installed for the following work.


1. Run the SIXABCon configuration program, see SIDEXIS installation of single/multistation systems. (Otherwise the service routines S.03, S.04 and S.05 cannot be performed.)
2. Install the Pan/Ceph sensor floppy using the SIXABCon configuration program.
3. Set SIDEXIS to factory service.

ATTENTION

IMPORTANT NOTES

- *For the following tests, it is important to ensure compliance with the **radiation protection guidelines** (see Operating Instructions).*
- *The **radiation** remains on as long as the release button is pressed.*
- *As long as radiation is active the **X-RAY** radiation on/off indicator must be lit up.*
- *The radiation on-time is also indicated **by an acoustic signal** emitted by the unit.*
- *The prescribed **cooling periods** are monitored by the **automatic exposure blocking** (see Operating Instructions).*

33 (C) Generating a test image (Panoramic unit)

1. Reconnect the unit to the power supply and switch it on.

2. Wait 1 minute.

The system's self-adjustment function runs.

NOTE

Do not press the test button on the multitimer before generating the test image.

3. Start SIDEXIS.

4. Activate the **Service** menu.

5. Activate the **Constancy test** menu item.

6. In the **X-ray components** menu, activate the **New** menu window.

7. Enter the service password.
(see SIDEXIS Service Manual).

The **New X-ray component** window will appear.

8. Fill in all the input fields.

9. Press the Return key to continue. If you do not know all the information, enter **Unknown** or **?** .

10. Press the **OK** button.

11. Activate the **XC XP exposure performance** icon.

This displays the **Select test type** selection dialog:

- Service exposure
- Acceptance test
- Constancy test
- Cancel


12. If another dialog is displayed:
Cancel
and activate **XC XP**.

13. From the **Select test type** dialog, select
Service exposure.

The selection dialog **Select service exposure** appears.

14. From this dialog, select **Digital test pattern (1)**.

(C)

15. Press the release button.

This triggers a cycle **with no** release of radiation.

A dialog box then appears.

If **OK** appears on the display, the data path is in order.

16. Press the **OK** button.

The test image generated is displayed on the screen.


If **Error** appears on the display, please refer to the Service Manual.

17. Close the test image ([Ctrl]+[F4]).

18. Press the Return key **R**.

The X-ray tube unit moves to its start position.

34 (C) Generating a test image (Ceph unit)

Checking the position of the patient diaphragm:

1. Set the primary diaphragm **30** on the diaphragm wheel (locking button **D**).
2. Press the Return key **R**.
 The X-ray tube assembly moves into the position for cephalography.
 The patient diaphragm moves into the position for asymmetric exposures.
3. Set the primary diaphragm **40** on the diaphragm wheel (locking button **D**).
4. Press the Return key **R**.
 The patient diaphragm moves into the position for symmetric exposures.

- If necessary adjust the diaphragms according to the installation instructions.

Generating a test image:

- From the **Select service exposure** dialog activate the option **Digital test pattern (1)**.

5. Press the release button.
 This triggers a scan with no release of radiation.

A dialog box then appears.
If

OK

appears on the display, the data path is OK.

Press the **OK** button.
The test image generated is displayed on the screen.


If

Error

appears on the display, see the information in the Service Manual.

- Close the test image.

35 (C) Phantom exposure

1. Insert the needle phantom as far as it will go.

2. Set the **primary diaphragm 10** on the diaphragm wheel.

3. Press the Return key **R**.

The X-ray tube unit moves to its start position.

4. From the selection dialog box on the screen, select **Service exposure** and activate **Quality test exposure (4)**.


5. Release an exposure for a **complete cycle** of the rotating part.

6. Evaluate the image on the screen.
Measure the line spacings using **Measure length** in the **Analyze** menu (see SIDEXIS Operator's Manual, Functional Description, Image processing, Measure length).

- If the tolerance given in the diagram is exceeded, you must adjust the M2/M3 actuators (see Service Manual).
- Exit exposure ([Ctrl]+[F4]).

7. Take the needle phantom back out again.

36 (C) Final Work

1. Screw the lead cover back on again.

ATTENTION

*For reasons of electromagnetic compatibility, **all** screws must be turned down and tightened.*

NOTE

Check that the side cover (59 68 250 D3200) has a recess for the LAN connection. If this is not the case, you must saw out/file the recess according to the above sketch.

2. Attach the side cover again using the two screws.

3. Reattach the front cover and screw it on tightly using the four screws.

4. Affix the nameplate with the new software version to the right side of the apron.

We reserve the right to make any alterations which may be due to technical improvements.

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Sirona Dental Systems GmbH

Fabrikstraße 31
D-64625 Bensheim
Germany
www.sirona.de

in the USA:
Sirona Dental Systems LLC
4835 Sirona Drive, Suite 100
Charlotte, NC 28273
USA

in Canada:
Sirona Canada
3250 Ridgeway Drive - Unit 5
Mississauga, Ontario L5L 5Y6
Canada

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© Sirona Dental Systems GmbH 2000